



# TRANSFER SETTINGS

MiniMed™ Paradigm® Veo™ Insulin Pump and Continuous Glucose Monitoring (CGM) System to MiniMed™ 770G pump

**STEP 1:** Copy your current settings onto the Insulin Pump Settings Form. If for any reason you cannot copy the settings, please contact your healthcare professional.

**STEP 2:** Insert a AA battery into your new pump. You can use either a lithium, alkaline, or rechargeable battery.

**STEP 3:** Follow the instructions found in this guide to enter the settings that you have recorded on the Insulin Pump Settings Form into your new pump. For additional information on your pump, please consult the MiniMed™ Paradigm® Veo™ Insulin Pump and Continuous Glucose (CGM) System User Guide.

**Medtronic**

We recommend that you record all settings on this form to refer to as you enter them into the new pump.

**BASAL SETTINGS**

**STANDARD BASAL RATE:**

Main ▶ Basal ▶ Basal Review

- 1: 12:00AM \_\_\_\_\_ U/H
- 2: \_\_\_\_\_ U/H
- 3: \_\_\_\_\_ U/H
- 4: \_\_\_\_\_ U/H
- 5: \_\_\_\_\_ U/H
- 6: \_\_\_\_\_ U/H

**MAX BASAL RATE:**

Main ▶ Basal ▶ Basal Setup ▶ Max Basal Rate

\_\_\_\_\_ U/H

**BASAL PATTERNS:**

Main ▶ Basal ▶ Basal Setup ▶ Patterns

On  Off

**Basal Pattern A:**

Main ▶ Basal ▶ Basal Review ▶ Pattern A

- 1: 12:00AM \_\_\_\_\_ U/H
- 2: \_\_\_\_\_ U/H
- 3: \_\_\_\_\_ U/H
- 4: \_\_\_\_\_ U/H

**Basal Pattern B:**

Main ▶ Basal ▶ Basal Review ▶ Pattern B

- 1: 12:00AM \_\_\_\_\_ U/H
- 2: \_\_\_\_\_ U/H
- 3: \_\_\_\_\_ U/H
- 4: \_\_\_\_\_ U/H

**UTILITIES**

**ALARM**

Main ▶ Utilities ▶ Alarm ▶

**Alert Type:**  Beep Long  Beep Medium  
 Beep Short  Vibrate

**Auto Off:**  On  Off \_\_\_\_\_ Hrs

**Low Resv Warning:**  \_\_\_\_\_ Units  \_\_\_\_\_ Time

**BOLUS SETTINGS**

**BOLUS WIZARD SETUP:**

Main ▶ Bolus ▶ Bolus Setup ▶ Bolus Wizard Setup ▶ Review Settings

**Wizard:**  On  Off **Carb Units:**  Grams  Exch

**Carb Ratios:** 1: 12:00AM \_\_\_\_\_ grams/U  
 2: \_\_\_\_\_ grams/U  
 3: \_\_\_\_\_ grams/U

**Sensitivity:** 1: 12:00AM \_\_\_\_\_ mmol/L/U  
 2: \_\_\_\_\_ mmol/L/U  
 3: \_\_\_\_\_ mmol/L/U

**BG Target:** 1: 12:00AM \_\_\_\_\_ - \_\_\_\_\_ mmol/L  
 2: \_\_\_\_\_ - \_\_\_\_\_ mmol/L  
 3: \_\_\_\_\_ - \_\_\_\_\_ mmol/L

**Active Ins Time:** \_\_\_\_\_ Hrs

**MAX BOLUS:**

Main ▶ Bolus ▶ Bolus Setup ▶ Max Bolus

\_\_\_\_\_ units

**SCROLL RATE:**

Main ▶ Bolus ▶ Bolus Setup ▶ Scroll Rate

0.1 U  0.05 U  0.025 U

**DUAL/SQUARE BOLUS:**

Main ▶ Bolus ▶ Bolus Setup ▶ Dual Square Bolus

On  Off

**EASY BOLUS:**

Main ▶ Bolus ▶ Bolus Setup ▶ Easy Bolus

On/Set  Off \_\_\_\_\_ units

**BG REMINDER:**

Main ▶ Bolus ▶ Bolus Setup ▶ BG Reminder

On  Off

**MISSED BOLUS REMINDER:**

Main ▶ Bolus ▶ Bolus Setup ▶ Missed Bolus Reminder

On  Off

1: \_\_\_\_\_

2: \_\_\_\_\_

**CONTINUOUS GLUCOSE MONITORING SETTINGS**

**LOW SUSPEND ALARM (MINIMED™ PARADIGM® VEO™ INSULIN PUMP AND CONTINUOUS GLUCOSE (CGM) SYSTEM ONLY)**

On \_\_\_\_\_ mmol/L (2.2 - 6.1 mmol/L)  Off

**LOW GLUCOSE ALERT**

Start on _____ (date)		
	<b>Time</b>	<b>Units</b>
	12 a.m.	_____ mmol/L
<input type="checkbox"/> On	_____	_____ mmol/L
	_____	_____ mmol/L
<b>OR</b>		_____ mmol/L
		<input type="checkbox"/> Off

**LOW PREDICTIVE ALERT**

On \_\_\_\_\_ minutes (5–30 min)  Off

**FALL RATE ALERT**

On \_\_\_\_\_ mmol/minute (0.065 - 0.275)  Off

**LOW REPEAT**

\_\_\_\_\_ minutes (5 min–1 hour)

Low Repeat is turned on to the default of 20 minutes when any Low Alert setting is programmed.

**HIGH GLUCOSE ALERT**

Start on _____ (date)		
	<b>Time</b>	<b>Units</b>
	12 a.m.	_____ mmol/L
<input type="checkbox"/> On	_____	_____ mmol/L
	_____	_____ mmol/L
<b>OR</b>		_____ mmol/L
		<input type="checkbox"/> Off

**HIGH PREDICTIVE ALERT**

On \_\_\_\_\_ minutes (5–30 min)  Off

**RISE RATE ALERT**

On \_\_\_\_\_ mmol/L/minute (0.065 - 0.275)  Off

**HIGH REPEAT**

\_\_\_\_\_ minutes (5 min–3 hours)

Default is set at 1 hour. May need to increase to 2 hours to minimize alerts.

## STARTUP WIZARD

Once you've inserted the battery, the pump powers up and the **Startup Wizard** begins, it will guide you through setting language, time format, and current time and date.

1. Select **English**.
2. When the **Select Time Format** screen appears, select a **12 Hour** or a **24 Hour** time format.
3. When **Enter Time** screen appears, the hour will be flashing. Press **^** or **v** to the correct hour and press **O**.
4. The minutes will be flashing. Press **^** or **v** to the correct minutes and press **O**.
5. The **AM / PM** will be flashing. Press **^** or **v** if needed and press **O**. Select **Next**.
6. When the **Enter Date** screen appears, select **Year**. Press **^** or **v** to the correct year and press **O**. Repeat this step to set up **Month** and **Day**. Select **Next**.
7. The rewinding screen will briefly appear, followed by a confirmation message. Select **Ok**.



## SETTING MAX BASAL

**Menu > Options > Delivery Settings**

1. Select **Max Basal/Bolus** and press **O**.
2. Select **Max Basal** to enter the **U/hr** amount. Use **^** and **v** to set the **Max Basal** amount. Press **O**.
3. Select **Save**.



## SETTING MAX BOLUS

**Menu > Options > Delivery Settings**


1. Select **Max Basal/Bolus** and press **O**.
2. Select **Max Bolus** to enter the **U/hr** amount. Use **^** and **v** to set the **Max Bolus** amount. Press **O**.
3. Select **Save**.



## SETTING BASAL PATTERN

**Menu > Basal > Delivery Settings > Basal Pattern Setup**

1. Select **Basal 1**.
2. Select **Options**, then select **Edit**.
3. Press **O** to enter **End** time. Press **O** to leave the **End** time at **12:00A**. If you are setting multiple rates use **^** or **v** to change the **End** time and press **O**.
4. Press **^** to enter **U/hr** and press **O**.
5. Continue adding **End** times and basal rates as needed, otherwise select **Done**.
6. Review the basal rates. If you need to correct, press **⬅** to go back to edit mode. If Basal rates are correct select **Save**.













## TURN BOLUS WIZARD ON AND SETUP

Menu > Options > Delivery Settings > Bolus Estimate Setup

1. Select **Bolus Wizard** to turn **On**.
2. Press **✓** to see the rest of the message then select **Next**.
3. Review the description of **Carb Ratio** and select **Next**.
4. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple Carb Ratios use **▲** or **▼** to change the **End** time, then press **○**.
5. Press **▲** or **▼** to enter g/U then press **○**.  
Continue adding **End** times and **Carb Ratios** as needed.
6. Select **Next** when done.
7. Review the description of **Sensitivity Factor** and select **Next**.
8. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple **Sensitivity Factors**, use **▲** or **▼** to change the **End** time, then press **○**.
9. Press **▲** or **▼** to enter **mmol/L per U**, then press **○**.  
Continue adding **End** times and **Sensitivity Factors** as needed.
10. Select **Next** when done.
11. Review the description of **BG Target** and select **Next**.
12. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple **BG Targets**, use **▲** or **▼** to change the **End** time, then press **○**.
13. Press **▲** or **▼** to enter **Lo** and then press **○**.
14. Press **▲** or **▼** to enter **Hi** and then press **○**.  
Continue adding **End** times and **BG Targets** as needed.
15. Select **Next** when done.
16. Review the description of **Active Insulin Time** and select **Next**.
17. Press **○** to enter **Duration** time.
18. Press **▲** or **▼** to enter **hr** then press **○**.
19. Select **Save**.

The **Bolus Wizard** setup is now complete.




**Note:** The settings shown are for illustration purposes only — your settings may be different.

## TURN SENSOR ON

Menu > Options > Utilities > Sensor Settings

1. Select **Sensor** to turn feature **On**.  
You can now access the **SmartGuard** menu and enter the settings.



## PROGRAMMING HIGH SETTINGS

Menu > Options > SmartGuard

1. Select **High Setup**.
2. Press **○** on the time segment. Press **○** again to leave the **End** time at **12:00A**. If setting multiple time segments, press **▲** to enter **End** time of first segment, and press **○**.
3. Press **▲** or **▼** to set **Hi** limit and press **○**. Press **○** to continue.
4. Select each feature you wish to turn on. If a feature is on, select it again to turn it back off.
5. Once settings are selected, select **Next** and then select **Done**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press **▼** to **Snooze** and press **○**.  
Select **High Snooze**. Press **▲** or **▼** to enter the correct time and press **○**. Select **Save**.




Your **High Setup** is now complete.

## PROGRAMMING LOW SETTINGS

Menu > Options > SmartGuard

1. Select **Low Setup**.
2. Press **○** on the time segment. Press **○** again to leave the **End** time at **12:00A**. If setting multiple time segments, press **▲** to enter **End** time of first segment, and press **○**.
3. Press **▲** or **▼** to set **Lo** limit and press **○**. Press **○** to continue.
4. Select each feature you wish to turn on. If a feature is on, select it again to turn it back off.
5. Once settings are selected, select **Next** and then select **Done**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press **▼** to **Snooze** and press **○**.  
Select **Low Snooze**. Press **▲** or **▼** to the correct time and press **○**. Select **Save**.



Your **Low Setup** is now complete.


Turn to page 5 to continue programming sensor settings.

## TURN AUTO MODE ON

Menu > Options > SmartGuard

1. Select **Auto Mode**.
2. Select **Auto Mode** again to turn Auto Mode **On**.
3. Check the screen to make sure that Auto Mode is set to **On**.  
Select **Save**


**Note:** If SmartGuard Auto Mode is not yet ready, after selecting **Save**, an alert occurs instructing to check the Auto Mode Readiness screen.  
**0-284-4416**




## PAIR THE TRANSMITTER

Menu > Options > Utilities > Device Options > Pair Device

1. Attach your transmitter to the charger and place next to the pump.
2. Select **Search** on the pump and immediately remove the transmitter from the charger.
3. The **Select Device** screen appears with a list of available devices. Select the **CGM** device that matches the serial number (**SN**) on your transmitter.
4. Make sure the serial number on the pump screen and transmitter match and select **Confirm**.

Your pump displays a message if the pump and transmitter are paired successfully. If the Sensor feature is turned on, the Connection icon  appears on the Home screen.



Please take some time to review ALL your entries and ensure they are correct.

## In Need of Technical Assistance?

**24-Hour Technical Support** Is Here to Help

**Call 1-800-284-4416**

For medical advice, contact your healthcare team.

**Note:** If you turn Auto Mode on before the 48 hour warm up is complete, you may receive an alert at midnight to Enter a BG for Auto Mode. To avoid this, you may want to wait to turn on Auto Mode on day 3 of your pump wear.