Airport Information

Medtronic

Airport security
- The pump must not go through the x-ray machine that is used for carry-on or checked luggage. - The full body scanner is also a form of x-ray. Remove the insulin pump and disconnect the infusion set at the insertion site before going through the fully body scanner. If continuous glucose monitoring (CGM) is used, remove the sensor and transmitter before the scan. - To avoid removing the devices, request an alternative screening process that does not use x-ray. - The insulin pump, infusion set, reservoir, and CGM system can withstand exposure to airport metal detectors used at airport security checkpoints.
In flight

The MiniMed™ insulin pump system is suitable for use in aircraft. Always monitor your blood glucose (BG) during air travel. Changes in air pressure that occur during flight takeoff and landing can cause over-delivery or under-delivery of insulin, which may result in hypoglycemia or hyperglycemia. Be ready to respond to alerts and symptoms.

Patient information
Name:
Address:
Country:
Phone number (include country code):

Medical device information
Medical device type:
Device serial number:

© 2025 Medtronic and Medtronic logo are trademarks of Medtronic.
™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

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Medical Emergency Card

I wear a required medical device prescribed by my physician.

Medtronic

In an emergency, please contact
Physician name:
Phone number:
Emergency contact:
Phone number:

I have diabetes
If my behavior is peculiar, if I appear intoxicated, or if I am unconscious, it may be a result of severe low blood sugar.
I am not intoxicated Call for medical assistance

If I am awake and able to swallow, give me a source of sugar (for example: juice, candy, or non-diet soft drink). Do not try to give me food or drink if I am unconscious.

Healthcare professionals
I wear an insulin pump.
- The pump delivers rapid-acting insulin at a constant rate. - To stop insulin delivery on the pump: - Press any button to wake the pump. - From the home screen, press ^ to go to the status screen. - Press ☉ to Suspend All Delivery. - Press > to highlight Yes. - Press ☉ to suspend all insulin delivery. Do not remove the pump battery. Do not remove the pump without medical consent. - If the pump sounds an alarm, follow the instructions on the pump screen or contact a local Medtronic support representative for further assistance.