Simplera Sync

Sensor User Guide

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M063212C003_1

Contacts:

Africa:

Medtronic South Africa and Southern Africa
Office Reception Tel:
+27(0) 11 260 9300
Diabetes: 24/7 Helpline:
0800 633 7867
Sub-Sahara 24/7 Helpline:
+27(0) 11 260 9490

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O.E.S. Distrimed sh.p.k.
Tel.: +355 402270660

Argentina:

Corpomedica S.A.
Tel.: +(11) 4 814 1333
Medtronic Directo 24/7:
+0800 333 0752

Armenia (AM):

Exiol LLC
Tel.: +374 98 92 00 11
or +374 94 38 38 52

Australia:

Medtronic Australasia Pty. Ltd.
Tel.: 1800 668 670

Austria / Österreich:

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Medtronic Ltd.
Tel.: +44 1923-205167

USA:

Medtronic Diabetes Global Headquarters
24-Hour Technical Support: +1-800-646-4633
To order supplies: +1-800-843-6687

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Introduction

The Simplera Sync™ sensor (MMT-5120C4 and MMT-5120D4) with Bluetooth® wireless technology is a component of a Medtronic insulin pump system.

The Simplera Sync sensor converts small amounts of glucose from the interstitial fluid under the skin into an electronic signal. The sensor uses that signal to provide sensor glucose (SG) values to a compatible Medtronic insulin pump system.

Intended use

The Simplera Sync sensor is a single-patient, single-use component of a personal continuous glucose monitoring (CGM) system. It is intended to communicate via Bluetooth Low Energy (BLE) with a compatible Medtronic insulin pump system to provide glucose information for diabetes management. It calculates glucose concentrations based on collected signals from the interstitial fluid and transmits glucose and device data to the networked device.

The Simplera Sync sensor is intended for use in home and professional healthcare environments.

Indications for use

The Simplera Sync sensor is indicated for the management of diabetes. Refer to your compatible Medtronic insulin pump system for specific system indications.

The Simplera Sync sensor does not require calibration.

Contraindications

No contraindications are associated with Simplera Sync sensor use.

User safety

Warnings and precautions

Read this entire user guide before attempting to insert the Simplera Sync sensor. The inserter portion of the sensor does not work the same way as other Medtronic insertion devices. The sensor is not inserted the same way as other Medtronic sensors. Failure to follow directions may result in improper insertion, pain, or injury.

Do not use the Simplera Sync sensor adjacent to other electrical equipment that may cause interference with normal system operation. For more information on electrical equipment that may cause interference with normal system operation, see *Exposure to magnetic fields and radiation*, page 13.

Do not use continuous glucose monitoring if hydroxyurea, also known as hydroxycarbamide, is taken. Hydroxyurea is used to treat certain diseases, such as cancer and sickle cell anemia. Hydroxyurea use results in higher sensor glucose readings compared to blood glucose readings.

Taking hydroxyurea while using continuous glucose monitoring can result in inaccurate or missed alerts, and substantially higher sensor glucose readings in reports than actual blood glucose readings. Always check the label of any medication being taken to confirm if hydroxyurea or hydroxycarbamide is an active ingredient. If hydroxyurea is taken, consult a healthcare professional. Use additional blood glucose meter readings to verify glucose levels.

Always consult a healthcare professional before using sensor glucose values to make treatment decisions if a medication that contains acetaminophen or paracetamol is taken while wearing the sensor. Medications that contain acetaminophen or paracetamol can falsely raise sensor glucose readings. The level of inaccuracy depends on the amount of acetaminophen or paracetamol active in the body and can differ for each person. Falsely elevated sensor readings can result in over-administration of insulin, which can cause hypoglycemia. Medications that contain acetaminophen or paracetamol include, but are not limited to, cold medicines and fever reducers. Check the label of any medications being taken to see if acetaminophen or paracetamol is an active ingredient. Use additional blood glucose meter readings to confirm blood glucose levels.

Always examine the Simplera Sync sensor box for damage. If the sensor box is open or damaged, examine the sensor for damage. If the sensor is visibly damaged, discard the device to avoid possible contamination.

Do not use the Simplera Sync sensor if any part of the device is damaged. If the device is damaged, discard the device to avoid possible contamination.

Do not use the Simplera Sync sensor if the tamper band is broken, damaged, or missing from the device. The sensor is sterile and non-pyrogenic unless the device is damaged. If the tamper band is broken, damaged, or missing from the device, the sensor and the needle can be exposed to contamination. A sensor and needle exposed to contamination can cause site infection if inserted into the body.

Do not use the Simplera Sync sensor if the cap label is broken, damaged, or missing from the device. The sensor is sterile and non-pyrogenic unless the device is damaged. If the cap label is broken, damaged, or missing from the device, the sensor and the needle can be exposed to contamination. A sensor and needle exposed to contamination can cause insertion site infection if inserted into the body.

Do not unscrew or remove the Simplera Sync sensor cap until the device is ready to be used. Do not remove the cap and store the device for future use. The sensor is sterile and non-pyrogenic unless the cap is removed from the device or the tamper band is broken. If the cap is not on the device or the tamper band is broken, the sensor and the needle can be exposed to contamination. A sensor and needle exposed to contamination can cause site infection if inserted into the body.

Do not place the cap back onto the device after cap removal. Placing the cap back on the device could cause damage to the needle, prevent a successful insertion, and cause injury.

Do not change or modify the Simplera Sync sensor. Changing or modifying the sensor can result in improper insertion, pain, or injury.

Do not let children hold the Simplera Sync sensor without adult supervision. Do not let children put any part of the sensor in their mouth. This product poses a choking hazard for young children that can result in serious injury or death.

The safety of the Simplera Sync sensor system has not been studied in persons who are critically ill. The impact of medications common to these conditions on system performance is unknown.

The performance of the sensor has not been evaluated for persons on dialysis and may be unsafe under those conditions.

Watch for bleeding at the insertion site on top of the Simplera Sync sensor. If bleeding occurs, apply steady pressure with a sterile gauze pad or clean cloth placed on top of the sensor for up to three minutes. If bleeding continues, is significantly visible on top of the sensor, or if there is excessive pain or discomfort after insertion, follow these steps:

1. Remove the Simplera Sync sensor and continue to apply steady pressure until the bleeding stops.
2. Dispose of the Simplera Sync sensor. See *Disposal*, page 25.
3. Check the insertion site for redness, bleeding, irritation, pain, tenderness, or inflammation. If there is redness, bleeding, irritation, pain, tenderness, or inflammation, contact a healthcare professional.
4. Insert a new Simplera Sync sensor in a different location.

Some skin care products, such as sunscreens and insect repellents, can damage the Simplera Sync sensor. Do not allow skin care products to touch the sensor. Wash hands after using skin care products before touching the sensor. If any skin care products touch the sensor, immediately wipe the sensor with a clean cloth.

Serious incident reporting

If a serious incident related to the device occurs, immediately report the incident to Medtronic and to Health Canada.

Exposure to magnetic fields and radiation

Do not expose the Simplera Sync sensor to Magnetic Resonance Imaging (MRI) equipment, diathermy devices, or other devices that generate strong magnetic fields (for example, CT scan, or other types of radiation). Exposure to strong magnetic fields can cause the sensor to malfunction or be unsafe, or can result in serious injuries.

IEC 60601-1-2; Special EMC Precautions for Medical Electrical Equipment

1. Special Precautions regarding Electromagnetic Compatibility (EMC): This body worn device is intended to be operated within a reasonable residential, domestic, public or work environment where common levels of radiated "E" (V/m) or "H" fields (A/m) exist, such as cellular phones, Wi-Fi™, Bluetooth wireless technology, electric can openers, microwave and induction ovens. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the provided instructions, may cause harmful interference to radio communications.
2. Portable and mobile RF communications equipment can affect medical electrical equipment. If you encounter RF interference from a mobile or stationary RF transmitter, move away from the RF transmitter that is causing the interference.
3. Be careful when using the Simplera Sync sensor closer than 30 cm (12 in) to portable radio frequency (RF) equipment or electrical equipment. If the sensor must be used next to portable RF equipment or electrical equipment, observe the sensor to verify correct system operation. Degradation of the performance of the sensor could result.

Essential performance

The Simplera Sync sensor performs measurements and estimates blood glucose concentration, and wirelessly communicates glucose information to the compatible Medtronic insulin pump system.

Risks

General risks with Simplera Sync sensor use include the following:

- Skin irritation or other reactions
- Bruising
- Discomfort
- Redness
- Bleeding
- Pain
- Rash
- Infection
- Raised bump
- Appearance of a small freckle-like dot where needle was inserted
- Allergic reaction
- Fainting secondary to anxiety or fear of needle insertion
- Soreness or tenderness
- Swelling at insertion site

- Sensor filament fracture, breakage, or damage
- Minimal blood splatter associated with sensor needle removal
- Residual redness associated with adhesive or tapes or both
- Scarring

Allergens


The Simplera Sync sensor contains nickel in stainless steel.

Reagents

The Simplera Sync sensor contains two biological reagents: glucose oxidase, and human serum albumin (HSA).

Glucose oxidase is derived from **Aspergillus niger** and manufactured to meet industry requirements for extraction and purification of enzymes for use in diagnostic, immunodiagnostic, and biotechnical applications. The HSA used in the Simplera Sync sensor consists of purified and dried albumin fraction V, derived from pasteurized human serum which is cross-linked via glutaraldehyde. Approximately 3 µg of glucose oxidase and approximately 10 µg of HSA are used to manufacture each sensor.


Simplera Sync sensor device components

Where to insert the Simplera Sync sensor

The image that follows shows insertion sites for ages 7 years and older. Target the shaded areas shown in the image, and make sure that the insertion site has a sufficient amount of fat.

Ages 7 years and older

Back of upper arm

Insertion on abdomen and upper buttocks is not approved.

Note: Assistance from another person may be needed for sensor insertion into the back of the upper arm. If assistance is not needed, a mirror may be helpful for self-insertion.

Inserting the Simpler Sync sensor

Preparing for insertion

1


The inserter label is on the top of the inserter. Before insertion, perform the following steps:

- Inspect the expiration date. Do not use an expired Simpler Sync sensor.
- Make note of the serial number (SN) and the CODE. Both numbers will be used later to pair the sensor with a compatible insulin pump.

Note: The SN and CODE label is also on the inside of the Simpler Sync sensor box lid.


2

Inspect the cap label for damage before insertion.

Note: Do not use the Simplera Sync sensor if the cap label is damaged or missing from the cap.

3

Inspect the tamper band to make sure that it is not broken, damaged, or missing from the device.

Note: Do not use the Simplera Sync sensor if the tamper band is broken, damaged, or missing.


4

Wash hands thoroughly with soap and water.

Note: Wear gloves when inserting the Simplera Sync sensor into another person to avoid accidental contact with patient blood. Minimal bleeding can occur.

5


Choose an insertion site that has a sufficient amount of fat. For insertion sites see *Where to insert the Simplera Sync sensor*, page 16.

For the best sensor performance, and to avoid accidental sensor removal, do not insert the Simplera Sync sensor into the following areas:


- muscle, tough skin, or scar tissue
- areas that are constrained by clothing or accessories
- areas subjected to rigorous movement during exercise

6


Clean the insertion site with alcohol. Allow the insertion site to air dry.

7


Unscrew the cap from the inserter, breaking the tamper band.

Note: Do not use the Simplera Sync sensor if the tamper band is broken, damaged, or missing from the device.


Insertion

8

Place the inserter on top of the prepared insertion site.


9

Press the inserter firmly against the body until there is a click.


After insertion

10

Gently pull the inserter straight from the body.


11

Smooth down the sensor adhesive with a finger to ensure the sensor stays on the body for the entire duration of use.

Note: Use over-the-counter tape if desired for additional adhesion.

12

Pair the Simpler Sync sensor with a compatible insulin pump.

Note: The SN and CODE are required to pair the sensor with a compatible display device. For details on how to pair the sensor with a compatible insulin pump, see the compatible Medtronic insulin pump system user guide.

Bathing and swimming

While on the body, the sensor is protected against continuous immersion in water at a depth of 2.4 m (8 feet) for up to 30 minutes. Shower and swim without removing the sensor.

Removing the Simplera Sync sensor

To remove the Simplera Sync sensor:

1. Gently peel the sensor adhesive away from the body.
2. Dispose of the Simplera Sync sensor in accordance with all local laws and regulations. For additional information, see *Disposal*, page 25.

Cybersecurity

The Simplera Sync sensor is designed with cybersecurity protections to keep data secure on the device and when communicating with a compatible insulin pump or mobile app. These protections include: encryption, integrity checks, and authentication checks.

The Simplera Sync sensor uses Bluetooth[®] Low Energy (BLE) communication to connect to a compatible pump.

To keep the Simplera Sync sensor and data secure, follow the security guidelines provided.

- Do not share the CODE with others.
- Perform pairing with a compatible pump in a private location.
- Never use a sensor that has a broken case. Use of a sensor with a broken case may reduce cybersecurity protections for your therapy and data.

For questions or concerns about the cybersecurity of the Simplera Sync sensor refer to the Assistance section within this user guide or visit www.minimed.com/security. Users will be informed once the product has reached the end of life and end of support.

Simplera Sync sensor wireless communication

Quality of service

The Simplera Sync sensor connects to a compatible display device via a Bluetooth low-energy technology network. The sensor sends glucose data and system-related alerts to the compatible display device, which verifies the integrity of received data after wireless transmission. The quality of the connection is in accordance with the Bluetooth Specification v4.2.

Traveling by air

The Simplera Sync sensor is safe for use on commercial airlines. As travel rules are subject to change, it is advisable to check with your local airport security authority before traveling.

Guidance and manufacturer's declaration

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
Emissions Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	CISPR 11 Group 1, Class B	The transmitter uses RF energy only for system communications. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.


Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2 Test Level	IEC 60601-1-2 Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2, ±4, ±8, ±15 kV air	±8 kV contact ±2, ±4, ±8, ±15 kV air	For use in a typical domestic, commercial, or hospital environment.
Power frequency magnetic field IEC 61000-4-8	30 A/m	30 A/m	For use in a typical domestic, commercial, or hospital environment.
Proximity magnetic fields IEC 61000-4-39, Table 11	IEC 60601-1-2, Table 11	IEC 60601-1-2, Table 11	For use in a typical domestic, commercial, or hospital environment.


Guidance and Manufacturer's Declaration – Electromagnetic Immunity

Immunity Test	IEC 60601-1-2 Test Level	IEC 60601-1-2 Compliance Level	Electromagnetic Environment Guidance
Proximity fields from RF wireless communications equipment	IEC 60601-1-2, Table 9	IEC 60601-1-2, Table 9	For use in a typical domestic, commercial, or hospital environment. Portable and mobile RF communications equipment should be used no closer to any part of the transmitter than the recommended separation distance of 30 cm (12 in). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF electromagnetic fields IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.

Radiated power

	Effective radiated power (ERP)	0.36 mW (-4.45 dBm)
	Effective isotropic radiated power (EIRP)	0.59 mW (-2.30 dBm)

	Effective radiated power (ERP)	0.67 mW (-1.75 dBm)
	Effective isotropic radiated power (EIRP)	1.10 mW (0.40 dBm)

Maintenance

Operation

Operating temperature range	2 °C to 40 °C (36 °F to 104 °F)
Air pressure range	70.33 kPa to 106.17 kPa (10.2 psi to 15.4 psi)
Operating relative humidity (RH) range	15% to 95%

Specifications

Biocompatibility	Sensor: Complies with EN ISO 10993-1
Applied parts	Sensor; Type BF applied part
Sensor glucose assay method	Enzymatic amperometric electrochemical sensor
Sensor glucose measurement range	2.8 to 22.2 mmol/L (50 to 400 mg/dL)
Sensor power source	2 silver oxide batteries
Sensor data	At least 2 days
Sensor memory	At least 2 days (glucose values stored every 5 minutes)
Sensor operating frequency	2.402 to 2.480 GHz BLE
Modulation	Gaussian Frequency-Shift Keying (GFSK)
Sensor transmission range	Minimum of 6.09 m (20 feet) line of sight in free-air
Sensor dimensions	2.865 × 2.865 × 0.477 cm (1.128 × 1.128 × 0.188 in)

Sensor weight

4.6 g (0.16 ounces)

Storage

CAUTION: Do not freeze the Simplera Sync sensor, or store it in direct sunlight, extreme temperatures, or high humidity. These conditions may damage the device.

Room temperature storage range	2 °C to 30 °C (36 °F to 86 °F)
Relative humidity (RH) storage range	Up to 95% relative humidity

Simplera Sync sensor life of use

The Simplera Sync sensor can be used one time and has a maximum life of up to 170 hours (7 days). The 170-hour life span of the sensor begins within 30 minutes after insertion.

CAUTION: Do not use the sensor if there is a sudden rise in sensor temperature. When operating the sensor in air temperatures of 40 °C (104 °F), under certain fault conditions, the temperature of the sensor may briefly rise up to 50 °C (121 °F). If there is a sudden rise in temperature or the sensor becomes hot or uncomfortable, remove and discard the sensor.

Disposal

Disposal requirements for electronic equipment, batteries, sharps and potential biohazardous materials differ based on location. Confirm disposal requirements for electronic equipment, batteries, sharps, and potential biohazardous materials with local laws and regulations.

The used inserter contains a needle which has been in contact with blood or other bodily fluids.

The used sensor contains a battery and has been in contact with blood or other bodily fluids. Disposal of the battery in any receptacle that could be exposed to extreme heat may cause the battery to ignite and result in serious injury.

Do not dispose of any component of this product with household waste or recyclables.

Dispose of the inserter and sensor in accordance with local laws and regulations.

Open Source Software (OSS) disclosure

This document identifies the Open Source Software that may be separately called, executed, linked, affiliated, or otherwise utilized by this product. Such Open Source Software is licensed to users subject to the terms and conditions of the separate software license agreement for such Open Source Software. Use of the Open Source Software shall be governed entirely by the terms

and conditions of such license. The source/object code and applicable license for the Open Source Software can be obtained at the following site: www.medtronicdiabetes.com/ossnotices.

Performance data

Refer to the Simplera Sync system user guide for performance data.

Assistance

Contact the local representative or refer to the local Medtronic website for assistance.

Refer to the Medtronic Diabetes International Contacts list in this user guide for contact information.

For definitions of the symbols displayed in the Simplera Sync sensor and package labels, see www.medtronicdiabetes.com/symbols-definitions.

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MiniMed and MiniMed logo are trademarks of Medtronic MiniMed, Inc. TM* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. The Bluetooth[®]* word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Medtronic is under license. Wi-FiTM* and Wi-Fi TM* logo are trademarks of Wi-Fi Alliance[®].

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© 2026 Medtronic